

Medical Device Alert

Immediate action update

Device

Silicone gel filled breast implants manufactured by Poly Implant Prothese (PIP).

All devices implanted after 01 January 2001.

Problem Action Implanting surgeons/implanting centres Clinical management of women implanted with • Identify women who were implanted with PIP silicone gel filled implants. PIP silicone gel filled implants after 01 January 2001. Reassure them that there is no current evidence of health risk associated with the Action by filler and there is no indication for routine action in the form of explantation or Medical directors. ultrasound. • Plastic surgeons and all surgeons involved in breast reconstruction and augmentation. Advise them that further information about • Directors of surgical units involved in breast the testing (see problem) is available on the MHRA website and that further reconstruction and augmentation · Nurse executive directors. information about clinical management is available on the websites of the British Specialist nurses involved in breast cancer Association of Plastic. Reconstructive and care. Aesthetic Surgeons (BAPRAS) and the Association of Breast Surgery (ABS). CAS deadlines **GPs** Action underway: 18 October 2010 Advise women who are concerned about Action complete: 29 November 2010 their PIP implants to consult their implanting surgeon or implanting centre. Note: the recommendations in this MDA update the advice given in MDA/2010/025.

Issued: 04 October 2010 at 16:30 Ref: MDA/2010/078

Problem

On Tuesday 30 March 2010 the French medical device regulatory authority (AFSSAPS) informed the MHRA that it had suspended the marketing, distribution, export and the use of silicone gel filled breast implants manufactured by PIP (a French breast implant manufacturer). AFSSAPS recalled all of these devices. The MHRA issued MDA/2010/025 on 31 March 2010 advising UK clinicians not to implant these devices.

AFSSAPS had carried out an inspection of the PIP manufacturing plant and established that breast implants manufactured by the company since 2001 had been filled with a silicone gel with a composition different from that approved.

AFSSAPS has carried out testing of affected implants to look at the genotoxicity (potential for cancer), chemical toxicity of the filler material and mechanical properties of the implant shell. The MHRA also commissioned tests to look at genotoxicity and chemical toxicity.

Test results have not shown any evidence of genotoxicity or chemical toxicity of the filler material. One of the French genotoxicity tests was, however, inconclusive and further testing will be conducted by AFSSAPS with results expected in early 2011.

Mechanical testing of the implant shells carried out by AFSSAPS suggests that there may be an increased risk of rupture. However, the Therapeutic Goods Administration (TGA) of the Australian government also carried out tests on PIP silicone gel filled breast implants and found that these implants conformed to the relevant international standards for this type of product, including those for gel cytotoxicity and shell strength.

Distribution

This MDA has been distributed to:

- NHS trusts in England (Chief Executives)
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Chief Executives)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

Onward distribution

Please bring this notice to the attention of all who need to know or be aware of it. This may include distribution by:

Trusts to:

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- Directors of surgical units involved in breast construction
- Medical directors
- Nurse executive directors
- Plastic surgeons and all surgeons involved in breast reconstruction or augmentation
- · Specialist nurses involved in breast cancer care

Primary care trusts to:

CAS liaison officers for onward distribution to all relevant staff including:

- General practitioners
- Practice managers
- Practice nurses

Establishments registered with the Care Quality Commission (CQC) (England only) to:

The MHRA considers this information to be important to:

- · Breast augmentation centres
- · Hospitals in the independent sector
- · Independent treatment centres
- · Private medical practitioners

Issued: 04 October 2010 at 16:30 Ref: **MDA/2010/078**

England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number MDA/2010/078 or 2010/003/030/081/019

Please note that telephone numbers for the MHRA contacts below will change on 25 October 2010. From that date please contact the MHRA Central Enquiry Point on 020 3080 6000 and ask for the person by name.

Technical aspects

Bayode Adisa or Ian Smith Tel: 020 7084 3223/3306

Email: bayode.adisa@mhra.gsi.gov.uk

ian.smith@mhra.gsi.gov.uk

Clinical aspects

Dr Susanne Ludgate Tel: 020 7084 3123

Email: susanne.ludgate@mhra.gsi.gov.uk

How to report adverse incidents

Please report via our website http://www.mhra.gov.uk

Further information about CAS can be found at https://www.cas.dh.gov.uk/Home.aspx

Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre

Health Estates Investment Group

Room 17 Annex 6

Castle Buildings

Stormont Estate

Dundonald BT4 3SQ

Tel: 02890 523 704

Fax: 02890 523 900

Email: NIAIC@dhsspsni.gov.uk

http://www.dhsspsni.gov.uk/index/hea/niaic.htm

How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website http://www.dhsspsni.gov.uk/niaic

Further information about SABS can be found at http://sabs.dhsspsni.gov.uk/

Issued: 04 October 2010 at 16:30 Ref: MDA/2010/078

Scotland

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre Health Facilities Scotland NHS National Services Scotland Gyle Square 1 South Gyle Crescent Edinburgh EH12 9EB

Tel: 0131 275 7575 Fax: 0131 314 0722 Email: nss.iric@nhs.net

http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-iric/

Wales

Enquiries in Wales should be addressed to:

Dr Sara Hayes Senior Medical Officer Medical Device Alerts Welsh Assembly Government Cathays Park Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: Haz-Aic@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health
© Crown Copyright 2010

Addressees may take copies for distribution within their own organisations