

# **Medical Device Alert**

# Device

Silicone tracheostomy tubes

Various manufacturers

Various models

| Problem   | Action  |
|---|---|
| The MHRA is aware of several incidents where the eyelets of tracheostomy tubes appear to have been damaged by the holder. This led to dislodgement of the tube with the subsequent risk of airway loss.  Action by  All staff responsible for the care of patients with tracheostomy tubes. | <ul> <li>Check the instructions for use for both the tracheostomy tube and the holder to ensure that they are compatible.</li> <li>Be aware that some manufacturers have recently updated their instructions for use to include a warning to users not to use their tracheostomy tubes with sharpedged holders, for example some types of Velcro or metal edged holders.</li> <li>Ensure that the instructions for use are followed for both the tracheostomy tube and the holder.</li> </ul> |
| CAS deadlines   |   |
| Action underway: 25 September 2012  |   |
| Action complete: 09 October 2012  |   |
| Note: These deadlines are for systems to be in place to take actions.   |   |

Issued: 11 September 2012 at 15:00 Ref: **MDA/2012/062** 

# **Device**

Tracheostomy tubes are held in position using holders (also known as neck ties). The holders are attached by being threaded through the eyelets on the flanges of the tracheostomy tube.

## Distribution

This MDA has been sent to:

- NHS trusts in England (Chief Executives)
- Care Quality Commission (CQC) (Headquarters) for information
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Equipment Co-ordinators)
- Local authorities in Scotland (Equipment Co-ordinators)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

#### **Onward distribution**

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

#### **Trusts**

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- · Adult intensive care units
- · All clinical departments
- All clinical staff
- All wards
- · Ambulance services directors
- · Ambulance staff
- Anaesthetists
- · Clinical governance leads
- · Day surgery units
- · Health and safety managers
- Hospital at home units
- Medical directors
- Medical libraries
- Nursing executive directors
- Outpatient departments
- Paediatric intensive care units
- Purchasing managers
- Risk managers
- Supplies managers
- Theatres

#### Primary care trusts

CAS liaison officers for onward distribution to all relevant staff including:

- · Community children's nurses
- · Community hospitals
- Community nurses
- District nurses
- · General practitioners
- Health visitors
- · Minor injury units
- NHS walk-in centres
- Practice nursesSchool nurses
- Walk-in centres

#### Social services

Liaison officers for onward distribution to all relevant staff including:

• In-house residential care homes

#### Independent distribution

#### Establishments registered with the Care Quality Commission (CQC) (England only)

This alert should be read by:

- Care homes providing nursing care
- · Care homes providing personal care
- Hospices
- · Hospitals in the independent sector
- Independent treatment centres
- · Private medical practitioners

Issued: 11 September 2012 at 15:00 Ref: MDA/2012/062

#### **Establishments registered with OFSTED**

This alert should be read by:

- · Educational establishments with beds for children
- · Residential special schools

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: safetyalerts@dh.gsi.gov.uk and requesting this facility.

# **England**

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number MDA/2012/062 or 2012/001/018/401/003

#### **Technical aspects**

Louise Mulroy and Gica Leclerc Medicines & Healthcare products Regulatory Agency Floor 4 151 Buckingham Palace Road London SW1W 9SZ

Tel: 020 3080 7344 and 6610

Fax: 020 8754 3965

Email: louise.mulroy@mhra.gsi.gov.uk

gica.leclerc@mhra.gsi.gov.uk

#### Clinical aspects

Jonathan Plumb
Medicines & Healthcare products Regulatory Agency
Floor 4
151 Buckingham Palace Road

Tel: 020 3080 6648 Fax: 020 8754 3965

London SW1W 9SZ

Email: jonathan.plumb@mhra.gsi.gov.uk

### How to report adverse incidents

Please report via our website http://www.mhra.gov.uk

Further information about CAS can be found at https://www.cas.dh.gov.uk/Home.aspx

# Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre

Health Estates Investment Group

Room 17 Annex 6 Castle Buildings Stormont Estate Dundonald BT4 3SQ

Tel: 02890 523 704 Fax: 02890 523 900

Email: NIAIC@dhsspsni.gov.uk

http://www.dhsspsni.gov.uk/index/hea/niaic.htm

Issued: 11 September 2012 at 15:00 Ref: MDA/2012/062

### How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website http://www.dhsspsni.gov.uk/niaic Further information about **SABS** can be found at http://sabs.dhsspsni.gov.uk/

# Scotland

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre Health Facilities Scotland NHS National Services Scotland Gyle Square 1 South Gyle Crescent Edinburgh EH12 9EB

Tel: 0131 275 7575 Fax: 0131 314 0722

Email: nss.iric@nhs.net

http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-iric/

# Wales

Enquiries in Wales should be addressed to: Improving Patient Safety Team Medical Directorate Welsh Government Cathays Park Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: Haz-Aic@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health
© Crown Copyright 2012

Addressees may take copies for distribution within their own organisations